



# SimplyGo, ready for a stroll

Quick start guide

**PHILIPS**

RESPIRONICS

# SimplyGo - Getting started



Select the mode you've been prescribed



Pulse mode SimplyGo detects each breath intake and delivers the oxygen delivery accordingly.



Continuous-flow mode Continuous oxygen delivery from 0.5 to 2 LPM.



Sleep mode Pulse mode with an increased detection sensitivity to sense breath intake more accurately at night. SimplyGo automatically switches to the last continuous flow setting used if no breath is detected after a period of time.

# Operating SimplyGo



1 **Charge** the battery or use the charger.



2 Insert SimplyGo into its **carrying case**.



3 Insert the **charged battery**.



4 **Connect** a standard single lumen **nasal cannula**. Make sure it is not pinched or kinked.  
Maximum length up to 9 m (30 ft).



5 **Turn on/off** SimplyGo by pressing the power button twice.



6 **Select the mode** by pressing the mode (**M**) button. Continue pressing the Mode button to scroll through the *Pulse*, *Continuous flow*, or *Sleep* modes.



7 **Set the flow rate** in accordance with your oxygen prescription using the plus (+) or minus (-) buttons.



8 **Make sure you have all the accessories you need**. Install SimplyGo on the cart.



9 **Position the nasal cannula on your face** and breathe normally through the nose.



# SimplyGo on the plane

The United States Federal Aviation Administration (FAA) has approved the SimplyGo for in-flight use on commercial airliners.

Although the FAA and other international aviation regulatory bodies allow the use of approved portable oxygen concentrators on commercial airlines, many domestic and international carriers have additional requirements including advance notification, and possession of a physician's letter and oxygen prescription.

If you are planning to fly with your SimplyGo, contact your airline well in advance of travel.

## SimplyGo battery specification

Battery autonomy		
	Pulse mode setting 2*	3 hr
	Continuous flow 2 LPM	54 min

\*Based on 20 breaths per minute

Charge time: 2-3 hours per battery when fully discharged



Ready. Set. [SimplyGo!](#)